

CNTN4 Rabbit mAb

Catalog No: #56991

Package Size: #56991-1 50ul #56991-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CNTN4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	CNTN4 Antibody detects endogenous levels of total CNTN4
Immunogen Description	A synthesized peptide derived from human CNTN4
Other Names	AXCAM; BIG 2; Cntn4; CNTN4A; contactin 4; SCA16;
Accession No.	Uniprot:Q8IWW2
Calculated MW	150kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of CNTN4 expression in Human cerebellum cell lysate.

Product Description

Contactins mediate cell surface interactions during nervous system development. Has some neurite outgrowth-promoting activity. May be involved in synaptogenesis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.