

CAP1 Rabbit mAb

Catalog No: #56995

Package Size: #56995-1 50ul #56995-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CAP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	CAP1 Antibody detects endogenous levels of total CAP1
Immunogen Description	A synthesized peptide derived from human CAP1
Other Names	CAP1; CAP1 PEN;
Accession No.	Uniprot:Q01518
Calculated MW	55kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500~1:2000

IHC: 1:50~1:200

ICC/IF: 1:50~1:200

Images



Western blot analysis of CAP1 expression in Raji cell lysate.

Product Description

Directly regulates filament dynamics and has been implicated in a number of complex developmental and morphological processes, including mRNA localization and the establishment of cell polarity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.