

PPP6C Rabbit mAb

Catalog No: #56998

Package Size: #56998-1 50ul #56998-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | PPP6C Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB ICC/IF FC |
| Species Reactivity | Human |
| Specificity | PPP6C Antibody detects endogenous levels of total PPP6C |
| Immunogen Description | A synthesized peptide derived from human PPP6C |
| Other Names | PP6; PP6C; PPP6; Ppp6c; |
| Accession No. | Uniprot:O00743 |
| Calculated MW | 32kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000IC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of PPP6C expression in HeLa cell lysate.

Product Description

Catalytic subunit of protein phosphatase 6 (PP6). PP6 is a component of a signaling pathway regulating cell cycle progression in response to IL2 receptor stimulation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.