

CLPTM1 Rabbit mAb

Catalog No: #57000

Package Size: #57000-1 50ul #57000-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CLPTM1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	CLPTM1 Antibody detects endogenous levels of total CLPTM1
Immunogen Description	A synthesized peptide derived from human CLPTM1
Other Names	HS9; N14;
Accession No.	Uniprot:O96005
Calculated MW	90kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500~1:2000

IHC: 1:50~1:200

ICC/IF: 1:50~1:200

FC: 1:50

Images



Western blot analysis of CLPTM1 expression in HepG2 cell lysate.

Product Description

CLPTM1 may play a role in T-cell development. A chromosomal translocation involving CLPTM1 is found in a family with cleft lip and palate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.