

## DPF2 Rabbit mAb

Catalog No: #57013

Package Size: #57013-1 50ul #57013-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DPF2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	DPF2 Antibody detects endogenous levels of total DPF2
Immunogen Description	A synthesized peptide derived from human DPF2
Other Names	BAF45D; DPF2; REQ; UBID4;
Accession No.	Uniprot:Q92785
Calculated MW	44kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

## Images



Western blot analysis of DPF2 expression in HeLa cell lysate.

## Product Description

May be a transcription factor required for the apoptosis response following survival factor withdrawal from myeloid cells. Might also have a role in the development and maturation of lymphoid cells.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.