

TRPV5 Rabbit mAb

Catalog No: #57024

Package Size: #57024-1 50ul #57024-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TRPV5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	TRPV5 Antibody detects endogenous levels of total TRPV5
Immunogen Description	A synthesized peptide derived from human TRPV5
Other Names	CAT2; ECaC; ECaC1; OTRPC3; TrpV5;
Accession No.	Uniprot:Q9NQA5
Calculated MW	90kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of TRPV5 expression in HepG2 cell lysate.

Product Description

Constitutively active calcium selective cation channel thought to be involved in Ca²⁺ reabsorption in kidney and intestine.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.