

## SNAP29 Rabbit mAb

Catalog No: #57025

Package Size: #57025-1 50ul #57025-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SNAP29 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	SNAP29 Antibody detects endogenous levels of total SNAP29
Immunogen Description	A synthesized peptide derived from human SNAP29
Other Names	CEDNIK; SNAP29;
Accession No.	Uniprot:O95721
Calculated MW	27kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:50

## Images



Western blot analysis of SNAP29 expression in (1) HeLa cell lysate; (2) Mouse brain lysate.

## Product Description

Involved in multiple membrane trafficking steps.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.