

## CCDC98 Rabbit mAb

Catalog No: #57034

Package Size: #57034-1 50ul #57034-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	CCDC98 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human
Specificity	CCDC98 Antibody detects endogenous levels of total CCDC98
Immunogen Description	A synthesized peptide derived from human CCDC98
Other Names	ABRA1; ABRAXAS; FAM175A;
Accession No.	Uniprot:Q6UWZ7
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details


WB:1:500~1:2000

ICC/IF:1:50~1:200

IP:1:50

FC:1:50

## Images



Western blot analysis of CCDC98 expression in MCF7 cell lysate.

## Product Description

Component of the BRCA1-A complex, a complex that specifically recognizes 'Lys-63'-linked ubiquitinated histones H2A and H2AX at DNA lesions sites, leading to target the BRCA1-BARD1 heterodimer to sites of DNA damage at double-strand breaks (DSBs).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.