

## COPS4 Rabbit mAb

Catalog No: #57037

Package Size: #57037-1 50ul #57037-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |  |
|-----------------------|--|
| Product Name          | COPS4 Rabbit mAb   |
| Host Species          | Rabbit   |
| Clonality             | Monoclonal   |
| Isotype               | Rabbit IgG   |
| Purification          | Affinity-chromatography  |
| Applications          | WB   |
| Species Reactivity    | Human Mouse Rat  |
| Specificity           | COPS4 Antibody detects endogenous levels of total COPS4  |
| Immunogen Description | A synthesized peptide derived from human COPS4   |
| Other Names           | CH4; COPS4; CSN4; DCH4; SGN4;  |
| Accession No.         | Uniprot:Q9BT78   |
| Calculated MW         | 50kDa  |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of COPS4 expression in SH-SY5Y cell lysate.

## Product Description

Component of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.