

## HNF-4-alpha Rabbit mAb

Catalog No: #57045



Package Size: #57045-1 50ul #57045-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	HNF-4-alpha Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	HNF-4-alpha Antibody detects endogenous levels of total HNF-4-alpha
Immunogen Description	A synthesized peptide derived from human HNF-4-alpha
Other Names	HNF 4 alpha; HNF4; HNF4A; HNF4a7; HNF4a8; HNF4a9; Hnf4alpha; MODY; MODY1; NR2A1; NR2A21; TCF14; Tcf4;
Accession No.	Uniprot:P41235
Calculated MW	53kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of HNF-4-alpha expression in HepG2 cell lysate.

## Product Description

Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and intestine.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.