

## BLNK Rabbit mAb

Catalog No: #57064

Package Size: #57064-1 50ul #57064-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	BLNK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	BLNK Antibody detects endogenous levels of total BLNK
Immunogen Description	A synthesized peptide derived from human BLNK
Other Names	AGM4; BASH; Bca; BLNK; Ly57; Lyw57; SLP65;
Accession No.	Uniprot:Q8WV28
Calculated MW	70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200

## Images



Western blot analysis of BLNK expression in Raji cell lysate.

## Product Description

Functions as a central linker protein that bridges kinases associated with the B-cell receptor (BCR) with a multitude of signaling pathways, regulating biological outcomes of B-cell function and development. Plays a role in the activation of ERK/EPHB2, MAP kinase p38 and JNK.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.