

Wnt16 Rabbit mAb

Catalog No: #57069

Package Size: #57069-1 50ul #57069-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Wnt16 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Wnt16 Antibody detects endogenous levels of total Wnt16
Immunogen Description	A synthesized peptide derived from human Wnt16
Other Names	WNT16;
Accession No.	Uniprot:Q9UBV4
Calculated MW	41kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of Wnt16 expression in HepG2 cell lysate.

Product Description

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.