

## AKAP 95 Rabbit mAb

Catalog No: #57075

Package Size: #57075-1 50ul #57075-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AKAP 95 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	AKAP 95 Antibody detects endogenous levels of total AKAP 95
Immunogen Description	A synthesized peptide derived from human AKAP 95
Other Names	Akap8; AKAP95;
Accession No.	Uniprot:O43823
Calculated MW	95kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


## Application Details

WB: 1:500~1:2000

IHC: 1:50~1:200

ICC/IF: 1:50~1:200

## Images



Western blot analysis of AKAP 95 expression in HepG2 cell lysate.

## Product Description

Anchoring protein that mediates the subcellular compartmentation of cAMP-dependent protein kinase (PKA type II).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.