

## Pyrophosphatase 1 Rabbit mAb

Catalog No: #57094

Package Size: #57094-1 50ul #57094-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Pyrophosphatase 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	Pyrophosphatase 1 Antibody detects endogenous levels of total Pyrophosphatase 1
Immunogen Description	A synthesized peptide derived from human Pyrophosphatase 1
Other Names	inorganic pyrophosphatase 1; IOPPP; PP1; PPA1; PPase;Pyp;
Accession No.	Uniprot:Q15181
Calculated MW	33kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:50

## Images



Western blot analysis of Pyrophosphatase 1 expression in HepG2 cell lysate.

## Product Description

Belongs to the PPase family.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.