

## GPRC5B Rabbit mAb

Catalog No: #57095



Package Size: #57095-1 50ul #57095-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GPRC5B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	GPRC5B Antibody detects endogenous levels of total GPRC5B
Immunogen Description	A synthesized peptide derived from human GPRC5B
Other Names	Gprc5b; RAIG2;
Accession No.	Uniprot:Q9NZH0
Calculated MW	52kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of GPRC5B expression in Jurkat cell lysate.

## Product Description

Unknown. This retinoic acid-inducible G-protein coupled receptor provide evidence for a possible interaction between retinoid and G-protein signaling pathways.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.