

## TAGAP Rabbit mAb

Catalog No: #57125

Package Size: #57125-1 50ul #57125-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TAGAP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	TAGAP Antibody detects endogenous levels of total TAGAP
Immunogen Description	A synthesized peptide derived from human TAGAP
Other Names	ARHGAP47; IDDM21; Tagap; TAGAP1;
Accession No.	Uniprot:Q8N103
Calculated MW	65,90kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:50

## Images



Western blot analysis of TAGAP expression in Jurkat cell lysate.

## Product Description

May function as a GTPase-activating protein and may play important roles during T-cell activation.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.