

## PRPF4 Rabbit mAb

Catalog No: #57128

Package Size: #57128-1 50ul #57128-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PRPF4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	PRPF4 Antibody detects endogenous levels of total PRPF4
Immunogen Description	A synthesized peptide derived from human PRPF4
Other Names	hPrp4; HPRP4P; PRP4; Prp4p; PRPF4; RP70; SNRNP60;
Accession No.	Uniprot:O43172
Calculated MW	58kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

## Images



Western blot analysis of PRPF4 expression in HeLa cell lysate.

## Product Description

Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.