



Recombinant Chicken Avidin(AVD)

Catalog No: #AP89651

Package Size: #AP89651-1 20ug #AP89651-2 100ug #AP89651-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Chicken Avidin(AVD)
Brief Description	Recombinant Proteins
Host Species	Yeast
Purification	Greater than 85% as determined by SDS-PAGE.
Immunogen Description	Expression Region: 25-152aa Sequence Info:Full Length of Mature Protein
Target Name	His-tagged
Accession No.	P02701
Target Species	Gallus gallus (Chicken)
Calculated MW	15.8 kDa
Tag Info	C-terminal 6xHis-tagged
Target Sequence	ARKCSLTGKWTNDLGSNMTIGAVNSRGEFTGTYITAVTATSNEIKESPLHGTQNTINKRTQPTFGFTVNWKFS ESTTVFTGQCfidrngkevlktmwllrSSVNDIGDDWKATRVGINIFTRLRTQKE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.