

Chicken Ovomuroid

Catalog No: #AP89659



Package Size: #AP89659-1 1mg #AP89659-2 10mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chicken Ovomuroid
Brief Description	Recombinant Protein
Purification	>95% (SDS-PAGE)
Species Reactivity	Purified from Chicken egg white
Immunogen Description	Sequence Info:Full length protein
Other Names	Allergen Gal d I
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.