

# Chicken Ovotransferrin

Catalog No: #AP89668



Package Size: #AP89668-1 1mg #AP89668-2 10mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken Ovotransferrin
Brief Description	Recombinant Protein
Purification	>95% (SDS-PAGE)
Species Reactivity	Purified from Chicken egg white
Immunogen Description	Sequence Info:Full length protein
Other Names	Chicken Ovotransferrin <sup>o</sup> Allergen Gal d III <sup>o</sup> Conalbumin <sup>o</sup> Serum transferrin
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.