

Ovalbumin

Catalog No: #AP89680



Package Size: #AP89680-1 100mg #AP89680-2 1g

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Ovalbumin |
| Brief Description | Recombinant Protein |
| Purification | >95% (SDS-PAGE) |
| Immunogen Description | Sequence Info:Full length protein |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.