

Recombinant Human Protein kinase domain-containing protein, cytoplasmic(PKDCC), Escherichia Coli



Catalog No: #AP89652

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP89652-1 20ug #AP89652-2 100ug #AP89652-3 1mg

Description

Product Name	Recombinant Human Protein kinase domain-containing protein, cytoplasmic(PKDCC), Escherichia Coli
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Protein Length: Full length protein Expression Region: 1-493
Other Names	FLJ18197; Hypothetical protein BC007901; Hypothetical protein LOC91461; MGC125960; Pkdcc; PKDCC_HUMAN; Protein kinase domain containing cytoplasmic; Protein kinase domain containing, cytoplasmic homolog (mouse); Protein kinase domain-containing protein, cytoplasmic; Protein kinase-like protein SgK493; SGK493; Sugen kinase 493; Vertebrate lonesome kinase; VlK
Accession No.	Q504Y2
Tag Info	The tag type will be determined during production process
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.