

TTF-I Polyclonal Antibody

Catalog No: #33193



Package Size: #33193-1 50ul #33193-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TTF-I Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	TTF-I Polyclonal Antibody detects endogenous levels of TTF-I protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human TTF1.
Accession No.	Swiss Prot:Q15361GeneID:7270
Calculated MW	103kD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.