

## eIF3A Rabbit mAb

Catalog No: #57134

Package Size: #57134-1 50ul #57134-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	eIF3A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	eIF3A Antibody detects endogenous levels of total eIF3A
Immunogen Description	A synthesized peptide derived from human eIF3A
Other Names	EIF3; eIF3 theta; eIF3a; EIF3S10; P167; p180; p185; TIF32;
Accession No.	Q14152
Calculated MW	160kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of eIF3A expression in Jurkat cell lysate.

## Background

Component of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.