

CABP Rabbit mAb

Catalog No: #57140

Package Size: #57140-1 50ul #57140-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CABP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	CABP Antibody detects endogenous levels of total CABP
Immunogen Description	A synthesized peptide derived from human CABP
Other Names	CaBP1; Calbrain; Calcium binding protein 1; Calcium binding protein 5; Caldendrin; HCALB BR
Accession No.	Q9NZU7
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of CABP expression in U-87 MG cell lysate.

Background

Modulates calcium-dependent activity of inositol 1,4,5-triphosphate receptors (ITPRs). Inhibits agonist-induced intracellular calcium signaling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.