

Palladin Rabbit mAb

Catalog No: #57156

Package Size: #57156-1 50ul #57156-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Palladin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Palladin Antibody detects endogenous levels of total Palladin
Immunogen Description	A synthesized peptide derived from human Palladin
Other Names	CGI15; MYN; Palld; PNCA1; SIH002; Stx13bp1;
Accession No.	Q8WX93
Calculated MW	130kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Palladin expression in PC3 cell lysate.

Background

A synthesized peptide derived from human Palladin

Note: This product is for in vitro research use only and is not intended for use in humans or animals.