

# Somatostatin 28 Rabbit mAb

Catalog No: #57161



Package Size: #57161-1 50ul #57161-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Somatostatin 28 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Applications	IHC
Species Reactivity	Human Mouse Rat
Specificity	Somatostatin 28 Antibody detects endogenous levels of total Somatostatin 28
Other Names	Growth hormone release inhibiting factor; SMST; Somatostatin 28; Somatostatin; SST;
Accession No.	P61278
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

IHC 1:50~1:200

## Background

The cyclic tetradecapeptide somatostatin 28 is widely distributed throughout the body and is an important regulator of endocrine and nervous system function. It exerts its biologic actions by binding to specific high affinity receptors on the cell surface.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.