

# CD163L1 Rabbit mAb

Catalog No: #57167



Package Size: #57167-1 50ul #57167-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD163L1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Applications	IHC
Species Reactivity	Human
Specificity	CD163L1 Antibody detects endogenous levels of total CD163L1
Other Names	CD163b; CD163L1; CD164L1; M160; Scart1;
Accession No.	Q9NR16
Calculated MW	160kDa
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

IHC 1:50~1:200

## Background

Isoform 1 is highly expressed in the spleen, lymph nodes, thymus, and fetal liver and weakly expressed in bone marrow and no expression was found in peripheral blood leukocytes. Isoform 1 expression is restricted to the monocyte and macrophage cell lines. Isoform 2 is only expressed in spleen.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.