

Recombinant Human Secreted frizzled-related protein 5(SFRP5),Escherichia Coli

Catalog No: #AP89653

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89653-1 20ug #AP89653-2 100ug #AP89653-3 1mg

Description

Product Name	Recombinant Human Secreted frizzled-related protein 5(SFRP5),Escherichia Coli
Brief Description	Recombinant protein
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Immunogen Description	Expression Region: 30-317aa Sequence Info:Full Length of Mature Protein
Target Name	His-tagged
Other Names	SFRP5; FRP1B; SARP3; Secreted frizzled-related protein 5; sFRP-5; Frizzled-related protein 1b; FRP-1b; Secreted apoptosis-related protein 3; SARP-3
Accession No.	Q5T4F7
Target Species	Homo sapiens (Human)
Tag Info	N-terminal His-tagged
Target Sequence	E EYDYYGWQAE PLHGRSYSKP PQCLDIPADL PLCHTVGYKR MRLPNLLEHE SLAEVKQQAS SWLPLLAKRC HSDTQVFLCS LFAPVCLDRP IYPCRSLCEA VRAGCAPLME AYGFPWPEML HCHKFPLDND LCIAVQFGHL PATAPPVTKI CAQCEMEHSA DGLMEQMCSS DFVVKMRIKE IKIENGDRKL IGAQKKKKLL KPGPLKRKDT KRLVLHMKNG AGCPCQLDS LAGSFLVMGR KVDGQLLLMA VYRWDKKNKE MKFAVKFMFS YPCSLYYPFF YGAAEPH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.