

Histone H3 (Phospho-Thr3) Conjugated Antibody

Catalog No: #C14135



Package Size: #C14135-AF350 100ul #C14135-AF405 100ul #C14135-AF488 100ul #C14135-AF555 100ul #C14135-AF594 100ul #C14135-AF647 100ul #C14135-AF680 100ul #C14135-AF750 100ul #C14135-Biotin 100ul #C14135-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Histone H3 (Phospho-Thr3) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Phospho-Histone H3 (T3) Antibody detects endogenous levels of total Phospho-Histone H3 (T3)
Immunogen Description	A synthesized peptide derived from human Phospho-Histone H3 (T3)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	H3 histone family, member A; H3/A; H31; H3FA; H3FB; H3FC; H3FD; H3FF; H3FH; H3FI; H3FJ; H3FK; H3FL; HIST1H3A
Accession No.	Uniprot:P84243
Calculated MW	16kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.