

alpha Synuclein (Phospho-Ser129) Conjugated Antibody

Catalog No: #C14154

Orders: order@signalwayantibody.com

Package Size: #C14154-AF350 100ul #C14154-AF405 100ul #C14154-AF488 100ul #C14154-AF555 100ul #C14154-AF594 100ul #C14154-AF647 100ul #C14154-AF680 100ul #C14154-AF750 100ul #C14154-Biotin 100ul #C14154-Conjugated 50ul

Description

Product Name	alpha Synuclein (Phospho-Ser129) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Phospho-alpha Synuclein (S129) Antibody detects endogenous levels of total Phospho-alpha Synuclein (S129)
Immunogen Description	A synthesized peptide derived from human Phospho-alpha Synuclein (S129)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	NACP; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; SYN; SYUA; alpha-synuclein; PARK1; PARK4; PD1; SNCA;
Accession No.	Uniprot:P37840
Calculated MW	18kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.