

GATA3 (Phospho-Ser308) Conjugated Antibody

Catalog No: #C14167

Package Size: #C14167-AF350 100ul #C14167-AF405 100ul #C14167-AF488 100ul #C14167-AF555 100ul #C14167-AF594 100ul #C14167-AF647 100ul #C14167-AF680 100ul #C14167-AF750 100ul #C14167-Biotin 100ul #C14167-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GATA3 (Phospho-Ser308) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Phospho-GATA3 (S308) Antibody detects endogenous levels of total Phospho-GATA3 (S308)
Immunogen Description	A synthesized peptide derived from human Phospho-GATA3 (S308)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GATA-binding factor 3; Trans-acting T-cell specific transcription factor GATA-3; GATA3;
Accession No.	Uniprot:P23771
Calculated MW	48kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Transcriptional activator which binds to the enhancer of the T-cell receptor alpha and delta genes. Binds to the consensus sequence 5'-AGATAG-3'. Required for the T-helper 2 (Th2) differentiation process following immune and inflammatory responses.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.