

c-Myc (Phospho-Thr58) Conjugated Antibody

Catalog No: #C14184

Package Size: #C14184-AF350 100ul #C14184-AF405 100ul #C14184-AF488 100ul #C14184-AF555 100ul #C14184-AF594 100ul #C14184-AF647 100ul #C14184-AF680 100ul #C14184-AF750 100ul #C14184-Biotin 100ul #C14184-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	c-Myc (Phospho-Thr58) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Phospho-c-Myc (T58) Antibody detects endogenous levels of Phospho-c-Myc (T58)
Immunogen Description	A synthesized peptide derived from human Phospho-c-Myc (T58)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MRTL; MYC; Myc proto-oncogene protein; c-myc;
Accession No.	Uniprot:P01106
Calculated MW	57kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Myc a proto-oncogenic transcription factor that plays a role in cell proliferation, apoptosis and in the development of human tumors. Seems to activate the transcription of growth-related genes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.