

EGFR (Phospho-Tyr1068) Conjugated Antibody

Catalog No: #C14205

Package Size: #C14205-AF350 100ul #C14205-AF405 100ul #C14205-AF488 100ul #C14205-AF555 100ul #C14205-AF594 100ul #C14205-AF647 100ul #C14205-AF680 100ul #C14205-AF750 100ul #C14205-Biotin 100ul #C14205-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	EGFR (Phospho-Tyr1068) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Phospho-EGFR (Y1068) Antibody detects endogenous levels of Phospho-EGFR (Y1068)
Immunogen Description	A synthesized peptide derived from human EGFR
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1; ERBB; ERBB1; HER1; EGFR; Urogastrone;
Accession No.	Uniprot:P00533
Calculated MW	175kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

Receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses. Phosphorylation of EGFR at specific serine and threonine residues attenuates EGFR kinase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.