

c-Jun (Phospho-Thr91) Conjugated Antibody

Catalog No: #C14214



Package Size: #C14214-AF350 100ul #C14214-AF405 100ul #C14214-AF488 100ul #C14214-AF555 100ul #C14214-AF594 100ul #C14214-AF647 100ul #C14214-AF680 100ul #C14214-AF750 100ul #C14214-Biotin 100ul #C14214-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	c-Jun (Phospho-Thr91) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Phospho-c-Jun (T91) Antibody detects endogenous levels of Phospho-c-Jun (T91)
Immunogen Description	A synthesized peptide derived from human c-Jun
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	JUN; p39; V-jun avian sarcoma virus 17 oncogene homolog; Proto-oncogene c-Jun; AP1; Activator protein 1; Transcription factor AP-1;
Accession No.	Uniprot:P05412
Calculated MW	43kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

c-Jun is a member of the Jun family containing c-Jun, JunB, and JunD, and is a component of the transcription factor activator protein-1 (AP-1). AP-1 regulated genes exert diverse biological functions including cell proliferation, differentiation, and apoptosis, as well as transformation, invasion and metastasis, depending on cell type and context.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.