

B Raf (Phospho-Thr401) Conjugated Antibody

Catalog No: #C14231

Package Size: #C14231-AF350 100ul #C14231-AF405 100ul #C14231-AF488 100ul #C14231-AF555 100ul #C14231-AF594 100ul #C14231-AF647 100ul #C14231-AF680 100ul #C14231-AF750 100ul #C14231-Biotin 100ul #C14231-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	B Raf (Phospho-Thr401) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Phospho-B Raf (T401) Antibody detects endogenous levels of total Phospho-B Raf (T401)
Immunogen Description	A synthesized peptide derived from human Phospho-B Raf (T401)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	B-RAF; B-RAF proto-oncogene serine/threonine-protein kinase; BRAF1; RAFB1; RMIL; RMIL serine/threonine-protein kinase; c-RMIL; kinase B-Raf;
Accession No.	Uniprot:P15056
Calculated MW	85kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Product Description

Protein kinase involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron. Phosphorylates MAP2K1, and thereby contributes to the MAP kinase signal transduction pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.