

# Lyn (Phospho-Tyr396) Conjugated Antibody

Catalog No: #C14243

Package Size: #C14243-AF350 100ul #C14243-AF405 100ul #C14243-AF488 100ul #C14243-AF555 100ul #C14243-AF594 100ul #C14243-AF647 100ul #C14243-AF680 100ul #C14243-AF750 100ul #C14243-Biotin 100ul #C14243-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lyn (Phospho-Tyr396) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	Phospho-Lyn (Y396) Antibody detects endogenous levels of total Phospho-Lyn (Y396)
Immunogen Description	A synthesized peptide derived from human Phospho-Lyn (Y396)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	JTK8; LYN; Tyrosine-protein kinase Lyn;
Accession No.	Uniprot:P07948
Calculated MW	56kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Lyn a tyrosine kinase of the Src family that plays a role in B cell antigen-receptor signaling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.