

## Phosphotyrosine Conjugated Antibody

Catalog No: #C14247



Package Size: #C14247-AF350 100ul #C14247-AF405 100ul #C14247-AF488 100ul #C14247-AF555 100ul #C14247-AF594 100ul

#C14247-AF647 100ul #C14247-AF680 100ul #C14247-AF750 100ul #C14247-Biotin 100ul #C14247-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phosphotyrosine Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Phosphotyrosine Antibody detects endogenous levels of total Phosphotyrosine
Immunogen Description	A synthesized peptide derived from Phosphotyrosine
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Phospho-Tyrosine;
Accession No.	Uniprot:N/A
Calculated MW	kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Tyrosine phosphorylation plays a key role in cellular signaling. Research studies have shown that in cancer, unregulated tyrosine kinase activity can drive malignancy and tumor formation by generating inappropriate proliferation and survival signals. Antibodies specific for phospho-tyrosine have been invaluable reagents in these studies.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.