

Bad (Phospho-Ser112) Conjugated Antibody

Catalog No: #C14254

Package Size: #C14254-AF350 100ul #C14254-AF405 100ul #C14254-AF488 100ul #C14254-AF555 100ul #C14254-AF594 100ul #C14254-AF647 100ul #C14254-AF680 100ul #C14254-AF750 100ul #C14254-Biotin 100ul #C14254-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Bad (Phospho-Ser112) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Bad (S112) Antibody detects endogenous levels of total Phospho-Bad (S112)
Immunogen Description	A synthesized peptide derived from human Bad
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BAD; BBC2; BBC6; BCL2L8;
Accession No.	Uniprot:Q92934
Calculated MW	23kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Promotes cell death. Successfully competes for the binding to Bcl-X(L), Bcl-2 and Bcl-W, thereby affecting the level of heterodimerization of these proteins with BAX.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.