

# Lysophospholipase 1 Conjugated Antibody

Catalog No: #C14255

Package Size: #C14255-AF350 100ul #C14255-AF405 100ul #C14255-AF488 100ul #C14255-AF555 100ul #C14255-AF594 100ul #C14255-AF647 100ul #C14255-AF680 100ul #C14255-AF750 100ul #C14255-Biotin 100ul #C14255-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lysophospholipase 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Lysophospholipase 1 Antibody detects endogenous levels of total Lysophospholipase 1
Immunogen Description	A synthesized peptide derived from human Lysophospholipase 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	APT1; hAPT1; LPL1; LYPLA1; Lysophospholipase 1; LYSOPLA; LysoPLA I;
Accession No.	Uniprot:O75608
Calculated MW	25 kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Hydrolyzes fatty acids from S-acylated cysteine residues in proteins such as trimeric G alpha proteins or HRAS. Has depalmitoylating activity and also low lysophospholipase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.