

## FGFR3 (Phospho-Tyr724) Conjugated Antibody

Catalog No: #C14262



Package Size: #C14262-AF350 100ul #C14262-AF405 100ul #C14262-AF488 100ul #C14262-AF555 100ul #C14262-AF594 100ul #C14262-AF647 100ul #C14262-AF680 100ul #C14262-AF750 100ul #C14262-Biotin 100ul #C14262-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FGFR3 (Phospho-Tyr724) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	Phospho-FGFR3 (Y724) Antibody detects endogenous levels of total Phospho-FGFR3 (Y724)
Immunogen Description	A synthesized peptide derived from human Phospho-FGFR3 (Y724)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ACH; CD333; CEK2; FGFR3; Fibroblast growth factor receptor 3; HSGFR3EX; Hydroxyaryl protein kinase; JTK4; Tyrosine kinase JTK 4;
Accession No.	Uniprot:P22607
Calculated MW	96kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

## Product Description

Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.