

ErbB2 (Phospho-Tyr1139) Conjugated Antibody

Catalog No: #C14263

Package Size: #C14263-AF350 100ul #C14263-AF405 100ul #C14263-AF488 100ul #C14263-AF555 100ul #C14263-AF594 100ul #C14263-AF647 100ul #C14263-AF680 100ul #C14263-AF750 100ul #C14263-Biotin 100ul #C14263-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ErbB2 (Phospho-Tyr1139) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Phospho-ErbB2 (Y1139) Antibody detects endogenous levels of total Phospho-ErbB2 (Y1139)
Immunogen Description	A synthesized peptide derived from human Phospho-ErbB2 (Y1139)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD340; CerbB2; ERBB2; HER2; Herstatin; MLN19; NEU; NGL; p185erbB2; TKR1;
Accession No.	Uniprot:P04626
Calculated MW	180kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Protein tyrosine kinase that is part of several cell surface receptor complexes, but that apparently needs a coreceptor for ligand binding.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.