

# MCM2 (Phospho-Ser108) Conjugated Antibody

Catalog No: #C14276

Package Size: #C14276-AF350 100ul #C14276-AF405 100ul #C14276-AF488 100ul #C14276-AF555 100ul #C14276-AF594 100ul #C14276-AF647 100ul #C14276-AF680 100ul #C14276-AF750 100ul #C14276-Biotin 100ul #C14276-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MCM2 (Phospho-Ser108) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Phospho-MCM2 (S108) Antibody detects endogenous levels of total Phospho-MCM2 (S108)
Immunogen Description	A synthesized peptide derived from human Phospho-MCM2 (S108)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BM28; CCNL1; cdc19; CDCL1 Cyclin like 1; MCM2; Minichromosome maintenance protein 2;
Accession No.	Uniprot:P49736
Calculated MW	125kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Acts as component of the MCM2-7 complex (MCM complex) which is the putative replicative helicase essential for 'once per cell cycle' DNA replication initiation and elongation in eukaryotic cells.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.