

## Nucleolin (Phospho-Thr84) Conjugated Antibody

Catalog No: #C14279



Package Size: #C14279-AF350 100ul #C14279-AF405 100ul #C14279-AF488 100ul #C14279-AF555 100ul #C14279-AF594 100ul #C14279-AF647 100ul #C14279-AF680 100ul #C14279-AF750 100ul #C14279-Biotin 100ul #C14279-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nucleolin (Phospho-Thr84) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Phospho-Nucleolin (T84) Antibody detects endogenous levels of total Phospho-Nucleolin (T84)
Immunogen Description	A synthesized peptide derived from human Phospho-Nucleolin (T84)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	NCL; Nucl; Nucleolin; Protein C23;
Accession No.	Uniprot:P19338
Calculated MW	100kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Nucleolin is the major nucleolar protein of growing eukaryotic cells. It is found associated with intranucleolar chromatin and pre-ribosomal particles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.