

Nucleolin (Phospho-Thr76) Conjugated Antibody

Catalog No: #C14280



Package Size: #C14280-AF350 100ul #C14280-AF405 100ul #C14280-AF488 100ul #C14280-AF555 100ul #C14280-AF594 100ul #C14280-AF647 100ul #C14280-AF680 100ul #C14280-AF750 100ul #C14280-Biotin 100ul #C14280-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Nucleolin (Phospho-Thr76) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Phospho-Nucleolin (T76) Antibody detects endogenous levels of total Phospho-Nucleolin (T76)
Immunogen Description	A synthesized peptide derived from human Phospho-Nucleolin (T76)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MS1116; NCL; Nucl; Nucleolin; Protein C23;
Accession No.	Uniprot:P19338
Calculated MW	100kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Nucleolin is the major nucleolar protein of growing eukaryotic cells. It is found associated with intranucleolar chromatin and pre-ribosomal particles. It induces chromatin decondensation by binding to histone H1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.