

p130(Phospho-Thr986) Conjugated Antibody

Catalog No: #C14282



Package Size: #C14282-AF350 100ul #C14282-AF405 100ul #C14282-AF488 100ul #C14282-AF555 100ul #C14282-AF594 100ul #C14282-AF647 100ul #C14282-AF680 100ul #C14282-AF750 100ul #C14282-Biotin 100ul #C14282-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	p130(Phospho-Thr986) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse
Specificity	Phospho-p130(T986) Antibody detects endogenous levels of total Phospho-p130(T986)
Immunogen Description	A synthesized peptide derived from human Phospho-p130(T986)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	130 kDa retinoblastoma associated protein; IB4; p130; PRB2; Rb2; Rbl2; RBR2; Retinoblastoma like 2; Retinoblastoma-related protein 2; RIC128;
Accession No.	Uniprot:Q08999
Calculated MW	128kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Product Description

Key regulator of entry into cell division. Directly involved in heterochromatin formation by maintaining overall chromatin structure and, in particular, that of constitutive heterochromatin by stabilizing histone methylation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.