

JAK1 (Phospho-Tyr1034/Tyr1035) Conjugated Antibody

Catalog No: #C14286

Orders: order@signalwayantibody.com
 Package Size: #C14286-AF350 100ul #C14286-AF405 100ul #C14286-AF488 100ul #C14286-AF555 100ul #C14286-AF594 100ul

 Support: tech@signalwayantibody.com

#C14286-AF647 100ul #C14286-AF680 100ul #C14286-AF750 100ul #C14286-Biotin 100ul #C14286-Conjugated 50ul

Description

Product Name	JAK1 (Phospho-Tyr1034/Tyr1035) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Phospho-JAK1 (Y1034 + Y1035) Antibody detects endogenous levels of total Phospho-JAK1 (Y1034 + Y1035)
Immunogen Description	A synthesized peptide derived from human Phospho-JAK1 (Y1034 + Y1035)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	JAK 1; JAK 1A; JAK 1B; JAK1; JAK1A; JAK1B; JTK3;
Accession No.	Uniprot:P23458
Calculated MW	133kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Tyrosine kinase of the non-receptor type, involved in the IFN-alpha/beta/gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.