

# Flotillin 2 Conjugated Antibody

Catalog No: #C56002

Package Size: #C56002-AF350 100ul #C56002-AF405 100ul #C56002-AF488 100ul #C56002-AF555 100ul #C56002-AF594 100ul #C56002-AF647 100ul #C56002-AF680 100ul #C56002-AF750 100ul #C56002-Biotin 100ul #C56002-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Flotillin 2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Flotillin 2 Antibody detects endogenous levels of total Flotillin 2
Immunogen Description	A synthesized peptide derived from human Flotillin 2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	epidermal surface antigen; ESA; flotillin-2; Membrane component chromosome 17 surface marker 1; ESA1; M17S1
Accession No.	Uniprot:Q14254
Calculated MW	47kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles. May be involved in epidermal cell adhesion and epidermal structure and function.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.