

Pygopus 2 Conjugated Antibody

Catalog No: #C56011

Package Size: #C56011-AF350 100ul #C56011-AF405 100ul #C56011-AF488 100ul #C56011-AF555 100ul #C56011-AF594 100ul #C56011-AF647 100ul #C56011-AF680 100ul #C56011-AF750 100ul #C56011-Biotin 100ul #C56011-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Pygopus 2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse
Specificity	Pygopus 2 Antibody detects endogenous levels of total Pygopus 2
Immunogen Description	A synthesized peptide derived from human Pygopus 2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PP7910; PYGO2; Pygo2 protein; Pygopus 2;
Accession No.	Uniprot:Q9BRQ0
Calculated MW	41kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Pygopus 2 is involved in signal transduction through the Wnt pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.