

Beclin 1 Conjugated Antibody

Catalog No: #C56013

Package Size: #C56013-AF350 100ul #C56013-AF405 100ul #C56013-AF488 100ul #C56013-AF555 100ul #C56013-AF594 100ul #C56013-AF647 100ul #C56013-AF680 100ul #C56013-AF750 100ul #C56013-Biotin 100ul #C56013-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Beclin 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Beclin 1 Antibody detects endogenous levels of total Beclin 1
Immunogen Description	A synthesized peptide derived from human Beclin 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ATG6; ATG6 autophagy related 6 homolog; Bcl-2-interacting protein beclin; Beclin 1 autophagy related;Beclin-1;Beclin1; BECN1; Protein GT197; VPS30;
Accession No.	Uniprot:Q14457
Calculated MW	55KDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Plays a central role in autophagy. Required for the abscission step in cytokinesis. May play a role in antiviral host defense. Protects against infection by a neurovirulent strain of Sindbis virus.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.